

CLAIMS

What is claimed is:

5

A method of preconcentrating trace analytes by:
extracting polar and/or non-polar analytes through a solgel extraction medium.

1)10

2. A method according to claim 1 wherein said extracting step is further redefined as feeding a sample through a sol-gel coated inner surface of a tube and extracting the analytes from the sample with the sol-gel coating.

15

3. A method according to claim 2 wherein said feeding step is further defined as passing the sample through a capillary tube, the tube including a sol-gel coated inner surface.

4. A method according to claim 2 wherein said feeding step is further defined as passing the sample through a sol gel monolithic bed.

20

25

5. A method according to claim 1 wherein the organic component of the sol-gel is selected from the group including sol-gel-active forms and/or derivatives of poly(ethylene glycol), poly(methylphenylsiloxane), poly(dimethyldiphenylsiloxane), poly(dimethylgiloxane), octadecylsilane, octylsilane, crown ethers, cyclodextrins, calixarenes, dendrimers, poly(styrene), poly(styrene-divinylbenzene), poly(acrylate), molecularly imprinted polymers, etc.

30

6. A method according to claim 1 further including the step of thermally desorbing the analytes from the sol-gel extraction medium.

10

15

20

- 7. A method according to claim 1 further including the step of desorbing the analytes from the sol-gel extraction medium.
- 8. A method according to claim 6 further including the step of applying the extracted analytes to a GC capillary column.
 - 9. A method according to claim 7 further including the step of applying the extracted analytes to a liquid phase separation technique.
 - 10. A method according to claim 1 further including the steps of preconditioning the sol-gel prior to said extracting step.
 - 11. A method according to claim 8 wherein said preconditioning step is further defined as simultaneously heating and purging an inert gas over the sol-gel.
 - 12. A microextraction method including the steps of microextraction polar and non-polar analytes in a sol-gel extraction medium:
 - desorbing the analytes from the sol-gel and analyzing the extracted, desorbed analytes.